<Hospital_Name>

<Users_Name>

<Department>

<Customer_Address>

<Zip_Code> <City>

<Country_Name>

<Reference: 97422051-FA>

5 August 2025

Urgent Field Safety Notice - Urgent Medical Device Recall
Encore™ 26 Inflation Device
Encore™ 26 Advantage Kit
NephroMax™ Kit
UroMax Ultra™ Kit

Dear Materials Manager and/or Healthcare Provider,

Boston Scientific is initiating a removal of certain batches of the Encore™ 26 Inflation Device and certain batches of Encore™ 26 Advantage Kits, NephroMax™ Kits and UroMax Ultra™ Kits in which the impacted Encore 26 inflation devices are a component (see Appendix 1) due to the potential presence of foreign material.

The Encore 26 Inflation Device is used with balloon dilation catheters to create and monitor pressure in the balloon and to deflate the balloon. The recall is isolated to Encore 26 Inflation Device batches that contain a specific component which has the potential to produce foreign material. An internal investigation determined that during use, foreign material particles could migrate from affected Encore 26 Inflation Devices into a balloon dilatation catheter.

The Encore 26 Inflation Device may be used in a variety of clinical applications including interventional cardiology, neurovascular, endoscopy, urology and peripheral vascular procedures.

The most serious potential adverse health consequence associated with use in these different applications is as follows:

Cardiovascular	A life-threatening embolic event may occur, with a remote likelihood, through a cascade of events in which foreign material enters the balloon dilatation catheter and a balloon rupture occurs when already inserted in the patient.
Neurovascular	A life-threatening embolic event may occur, with a remote likelihood, through a cascade of events in which foreign material enters the balloon dilatation catheter and a balloon rupture occurs when already inserted in the patient.
Urology	Laceration/perforation requiring surgical intervention, with a remote likelihood, may result from deflation issues resulting from foreign material in the lumen of the balloon dilatation catheter.
Endoscopy	Obstruction at the biliary or pancreatic duct orifice leading to biliary obstruction or pancreatitis, inflammation, infection and additional intervention, with a remote likelihood, resulting from a cascade of events in which foreign material enters the balloon dilatation catheter and a balloon rupture occurs when already inserted in the patient.
Peripheral Vascular	An embolic event requiring surgical intervention may occur, with a remote likelihood, in which foreign material enters the balloon dilation catheter and a balloon rupture occurs when already inserted in the patient.

Boston Scientific has received 18 complaints of foreign material in the Encore device, a rate of 0.0009%, including one (1) complaint with the cascade of events of foreign material in the Encore device in conjunction with a balloon rupture, a rate of 0.000046%. No patient harm has been reported in any of these complaints.

In all clinical applications noted above, the most common potential adverse health consequence is prolongation of the procedure as the Encore inflation device is exchanged for another inflation device if the foreign material is detected during preparation or use of the device. There is a remote chance the foreign material could impact the ability to inflate or deflate the balloon used with the Encore 26 inflation device. Impact to balloon inflation may result in a prolongation of the procedure to exchange the device for another inflation device. Deflation issues resulting from foreign material in the lumen of the balloon may require the balloon to be ruptured so it can be removed from the patient. Harms that could result from removing a non-deflated balloon or needing to rupture the balloon for removal may include: infection, inflammation, obstruction, pain, laceration/perforation, hemorrhage, stenosis and unretrieved device fragments. There is a remote chance that any of these harms may occur. Boston Scientific has not received any reports of these patient harms.

Boston Scientific recommends that patients who have been treated with an Encore inflation device should continue to follow standard of care and any additional patient management should be at the discretion of their physician on an individualized basis.

Our records indicate that your facility has received some of the concerned product. **The table below (Appendix 1) provides a complete list of all affected products**, including Product Description, Material Number (UPN), GTIN and Lot/Batch numbers and expiry date. Please note that **only the devices listed below are affected. No other Boston Scientific product is involved in this Field Safety Notice. Further distribution or use of any remaining product affected by this action should cease immediately.**

INSTRUCTIONS:

1- **Please immediately discontinue use of the Boston Scientific product reported in the list and remove all of the affected units from your inventory**, regardless of where these units are stored in your facility. Segregate the units in a secure place, pending return to Boston Scientific.

2- **Please complete the attached Verification Form even if you do not have any product to return.**

3- **When completed, please return the Verification Form to your local Boston Scientific office** for the attention of <Customer_Service_Fax_Number> on or before **29 August 2025**.


4- **If you have products to return**, please package them in an appropriate shipping box. **After receipt of the Verification Form, Boston Scientific will contact you to arrange return. The entire Encore 26 Advantage Kit, NephroMax Kit, or UroMax Ultra Kit must be returned in order to be reimbursed.**

5- Please pass this notice to any health professional from your organization that needs to be aware and to any organization where the potentially affected devices have been transferred (if appropriate). Please provide Boston Scientific with details of any affected devices that have been transferred to other organizations (if appropriate).

Your Competent Authority is being notified of this Field Safety Notice.

Patient safety is Boston Scientific's highest priority. As such, we are committed to transparent communication with physicians and healthcare professionals to ensure you have timely, relevant information for managing your patients. If you require additional assistance or more information regarding this communication, please contact your local Boston Scientific representative.

Sincerely,

Kara Carter
Vice President, Global Quality

Enclosed: Appendix 1 Affected Product List
 Reply Verification Tracking Form

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